

EXTENSIN CABLE FOR IOP SENSOR (2M)

Ref. 554111



30/08/2024 20:50:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

EXTENSIN CABLE FOR IOP SENSOR (2M)



Ref. 554111

Features

CE Marking	CE	Medical device	No
Bar code	28809315676540	HS Code	9018 90 2090
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	2
Packed unit / Weight	0.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:50:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr