

REPLACEMENT HEAD SOFT GREEN FOR OTOSCOPE SMARTLED 5500 PEDIA



Ref. 570518P

 30/08/2024 20:15:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REPLACEMENT HEAD SOFT GREEN FOR OTOSCOPE SMARTLED 5500 PEDIA



Ref. 570518P

Features

CE Marking	NONE	Medical device	No
Bar code	3700446041755	HS Code	90189042
Country	Pakistan		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.046
Bare Unit / Weight (kg)	0.046		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:15:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr