

LED BULB FOR LITESTICK

Ref. 570690



30/08/2024 20:19:13 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LED BULB FOR LITESTICK



Ref. 570690

Description

Spare part for Litestick® penlight.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009908	HS Code	90189084
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:19:13 - For information purposes only, can be subject to change