

RECHARGING DOCK FOR SMARTLED 5500-R



Ref. 579534

 30/08/2024 20:22:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

RECHARGING DOCK FOR SMARTLED 5500-R



Ref. 579534

Features

CE Marking	NONE	Medical device	No
Bar code	3700446032913	HS Code	85366990
Country	Pakistan		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	8
Packed unit / Length (cm)	17	Packed unit / Weight	0.354
Bare Unit / Weight (kg)	0.354		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:22:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr