

CABLE USB-C - USB MALE 2M POUR SMARTLED 5500-R



Ref. 579535

 30/08/2024 20:22:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

**CABLE USB-C - USB MALE 2M POUR
SMARTLED 5500-R**

Ref. 579535

Features

CE Marking	NONE	Medical device	No
Bar code	3661040216782	HS Code	9018110000
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	8
Packed unit / Length (cm)	17	Packed unit / Weight	354
Bare Unit / Weight (kg)	354		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	80
Minimum humidity	5	Maximum humidity	90
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:22:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr