

PACK OF 4 STERILE GUEDEL CANNULAS SINGLE USE



Ref. 580840-



24/10/2024 20:45:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PACK OF 4 STERILE GUEDEL CANNULAS SINGLE USE



Ref. 580840-

Description

Set of 4 single-use cannulas in individual sterile sachet, one of each size: n°1/2/3/4
The Guedel's cannula is used to keep the airway of an unconscious patient open.

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446034665	HS Code	90183900
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.05
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

24/10/2024 20:45:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr