

DISPOSABLE TRANSDUCER (X50)

Ref. CSTRAN01

Sibelmed

CE
0318

DISPOSITIF MEDICAL
CLASSE
IIa
MEDICAL DEVICE



30/08/2024 23:07:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

DISPOSABLE TRANSDUCER (X50)



Ref. CSTRAN01

Description

Features

CE Marking	CE0318	Medical device	Class IIa
Bar code	3361040000045	HS Code	90181990
Country	Spain		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 23:07:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr