

## LARGE ADULT CUFF 32-43 CM FOR 4VS MONITOR



Ref. 551623



30/08/2024 20:01:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)

# LARGE ADULT CUFF 32-43 CM FOR 4VS MONITOR



Ref. 551623

## Description

Adult wide armband for PNI module of SPENGLER 4VS vital signs monitor.  
Dimensions (32-43cm)

## Features

|                             |               |                           |            |
|-----------------------------|---------------|---------------------------|------------|
| CE Marking                  | CE            | Medical device            | Class I    |
| Bar code                    | 3700446033071 | HS Code                   | 9018901000 |
| Warranty (months)           | 6             | Country                   | China      |
| DIMENSIONS AND WEIGHT       |               |                           |            |
| Packed unit / Width (cm)    | nc            | Packed unit / Height (cm) | nc         |
| Packed unit / Length (cm)   | nc            | Packed unit / Weight      | 0.088      |
| SPECIFIC STORAGE CONDITIONS |               |                           |            |
| Minimum temperature         | NA            | Max temperature           | NA         |
| Minimum humidity            | NA            | Maximum humidity          | NA         |
| OTHER FEATURES              |               |                           |            |
| Packed unit / Quantity      | 1             |                           |            |

 30/08/2024 20:01:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31  
 [info@gsh-med.fr](mailto:info@gsh-med.fr)