

NIPB TUBING 2,5M FOR 4VS MONITOR

Ref. 551624



 30/08/2024 20:01:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

NIPB TUBING 2,5M FOR 4VS MONITOR



Ref. 551624

Features

CE Marking	NONE	Medical device	No
Bar code	3661040213842	HS Code	85049099
Warranty (months)	6	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	16	Packed unit / Height (cm)	3
Packed unit / Length (cm)	17	Packed unit / Weight	0.118
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 20:01:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr