

MUCOSITY ASPIRATOR CA-MI NEW ASPIRET BOTTLE / BAG UU FLOVAC 1L



Ref. BCFLOV05



29/08/2024 21:17:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

MUCOSITY ASPIRATOR CA-MI NEW ASPIRET BOTTLE / BAG UU FLOVAC 1L



Ref. BCFLOV05

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

Description

Suitable for tracheostomy patients and is ideal for minor surgical applications as well as postoperative home treatments.

Main applications: ENT, Home Care, Bedroom, Tracheostomy

- Composed of: New Aspiret Suction Unit + Flovac Kit 1L (1x Flovac Jar 1L Reusable / 1x Sterile UU Pocket 1L Flovac/ 1x Set Tubes Ø6X10mm + connectors / 1x Round Retractor for Flovac System / 1x Yellow Male Connector FLOVAC System (5 pcs) / 1x Configuration Manual Suction Units)

- Optional

roller foot- The main system is equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Supplied accessories: Fluid collection jar with an anti-reflux valve with different options / 6 x 10 Silicone hoses (autoclavable) and conical connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko plug

- Engine: Oil-free piston pump requiring no maintenance

- Power supply: 230V-50Hz - other configurations available on request

- Maximum vacuum (adjustable): 0.75 Bar -75 kPa -563 mmHg

- Max free air flow: 15 l/min

- Noise level: 59, 6 dBA

- Operating cycle: Uninterrupted

- Weight: 2.2 kg / Dimensions: 23 x 19 x 16 cm

 29/08/2024 21:17:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MUCOSITY ASPIRATOR CA-MI NEW ASPIRET BOTTLE / BAG UU FLOVAC 1L



Ref. BCFLOV05

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040197562	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	23
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	16
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	19
Packed unit / Weight	2.2	Bare Unit / Weight (kg)	2.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 21:17:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr