

DIGITAL DERMATOSCOPE FIREFLY DE300A

Ref. DFIDE300A

Firefly

CE

DISPOSITIF MEDICAL
CLASSE
I
MEDICAL DEVICE

GARANTIE
1 AN
1 YEAR
GUARANTEE



30/08/2024 09:23:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

DIGITAL DERMATOSCOPE FIREFLY DE300A



Ref. DFIDE300A

Connected dermatoscope for telemedicine use

Description

- Resolution of 2 real megapixels without processing with interpolation techniques.
- Capture images or videos instantly: DE300 uses a state-of-the-art image processing engine that virtually eliminates visual lag, allowing for an intuitive user experience.
- Magnification up to 35x optically and 105x digitally
- Digital dermatoscope with built-in high quality polarizer offering multiple levels of polarization.
- Equipped with dual multilayer glass lenses, 650nm cut-off, thus offering superb depth of field and incomparable image sharpness.
- Designed for Digital Epiluminescence Microscopy (DELM)
- Easily integrates with EMR systems

Features

CE Marking	CE	Medical device	Class I
Bar code	8512890070718	HS Code	8525803000
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.292
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 09:23:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr