

SUNTECH CT50 BP STANDARD SPO2 PRINTER

SUNTECH

Ref. DSU222202

CE
0413

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
2 ANS
2 YEARS
GUARANTEE



 30/08/2024 09:25:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

SUNTECH CT50 BP STANDARD SPO2 PRINTER

SUNTECH

Ref. DSU222202

The world reference for the measurement of vital parameters. The CT 50 is a compendium of state-of-the-art technology

Description

The SUNTECH CT 50 monitor is a condensed version of the latest technology.

The world reference for measuring vital signs.

Nellcor SpO2 technology, optional HIS, USPx4, Ethernet LAN and WIFI connectivity.

3 blood pressure measurement modes: Automatic, Auscultatory and average mode.

8-inch colour touch screen, optional barcode reader, optional printer and optional wheeled base

Features

CE Marking	CE0413	Medical device	Class IIb
Bar code	3700446025120	HS Code	9018199560
Warranty (months)	24	Country	China

DIMENSIONS AND WEIGHT

Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	5.02

SPECIFIC STORAGE CONDITIONS

Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA

OTHER FEATURES

Packed unit / Quantity	1
------------------------	---

30/08/2024 09:25:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr