

MATRIX NEEDLE PRINTER FOR MWA

Ref. EKMWA06

KERN®



29/08/2024 17:30:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MATRIX NEEDLE PRINTER FOR MWA



Ref. EKMWA06

Description

- To print weighing values on normal paper, perfect accordingly for long-term archiving
- Needle printing
- Paper strips: 44 mm
- Case dimensions L××H 106×158×40 mm

Tip: Buy the scale and printer together

Features

CE Marking	NONE	Medical device	No
Bar code	3661040173696	HS Code	84231090
Warranty (months)	12	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	22	Bare Unit / Width (cm)	10.6
Packed unit / Height (cm)	13	Bare Unit / Height (cm)	4
Packed unit / Length (cm)	23	Bare Unit / Length (cm)	15.8
Packed unit / Weight	1	Bare Unit / Weight (kg)	0.48
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:30:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr