

WHEEL STAND FOR SPOT-CHECK 2VS

Ref. 546024



30/08/2024 19:54:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

WHEEL STAND FOR SPOT-CHECK 2VS



Ref. 546024

Description

Roller base for SPENGLER spot-check station 2VS

Features

CE Marking	NONE	Medical device	No
Bar code	3700446029715	HS Code	94029012003
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	10.8
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 19:54:12 - For information purposes only, can be subject to change