

CA-MI NEW ASKIR SUCTION UNIT 118 BOTTLE/BAG UU FLOVAC 2L



Ref. BCFLOV305



29/08/2024 21:06:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

CA-MI NEW ASKIR SUCTION UNIT 118 BOTTLE/BAG UU FLOVAC 2L



Ref. BCFLOV305

Electrical medical device for nasal and tracheal aspiration of body fluids in children or adults orally.

Description

CA-MI New Askir 118 Suction Unit - Bottle/Sachet UU Flovac® 2L

Electric medical device for nasal and tracheal aspiration of body fluids in children or adults orally.

Main applications: Emergency, Home Care, Patient Transfer, Room, Tracheostomy, Med. Veterinarian

Composed of: New Askir 118 Suction Unit + Flovac 2L Kit (1x Reusable Flovac 2L Flovac Jar / 1x Sterile UU Pocket 2L Flovac/ 1x Set tubes Ø6X10mm + connectors / 1x Round Retractor for Flovac® System / 1x Yellow Male Connector FLOVAC System (5 pcs) / 1x Configuration Manual Suction Units)

- Large LCD display for easy reading of vacuum values with vacuum adjustment keys to increase suction accuracy.

Supplied accessories: Fluid collection jar with anti-reflux valve with different options / 8x14 Silicone hoses (autoclavable) and conical connector 10-11-12 mm- Antibacterial & Hydrophobic filter (single patient) + CH20 cannula / Power cord, 12V automotive adapter, universal AC/DC adapter - 3 operating modes: mains 100-240V

, 12V

car adapter or and rechargeable internal battery make it a device that finds its Application in several areas and conditions.

- Long battery life and low noise level during operation.

PROXIMITY function is used to activate or deactivate the device without tactile contacts and prevents cross-contamination between two patients.

- Long-lasting brushless motor.

No smell of coal residues

Maximum vacuum (adjustable): -75 kPa -563 mmHg -0.75 Bar

Max free air flow: Hg 26 l/m (without jar)

)- Operating

ie: Uninterrupted (connected to the mains)

29/08/2024 21:06:26 - For information purposes only, can be subject to change

Weight: 2.7 kg / Dimensions: 25 x 19 x 15 cm

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CA-MI NEW ASKIR SUCTION UNIT 118 BOTTLE/BAG UU FLOVAC 2L



Ref. BCFLOV305

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040197241	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	2.7
Bare Unit / Weight (kg)	2.7		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			

29/08/2024 21:06:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr