

STETHOSCOPE II SPENGLER DISPLAY

Ref. K90026



30/08/2024 23:53:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

STETHOSCOPE II SPENGLER DISPLAY



Ref. K90026

Dare the design display to highlight your stethoscopes.

Description

Modular display that adapts to stethoscopes: Magister II, Cardioprestige II and Pulse II.

Features

CE Marking	NONE	Medical device	No
Country	France		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 23:53:59 - For information purposes only, can be subject to change