

EMERGENCY BAG LIGHT (22L / 44X25X27)

EMERGENCY'S™

Ref. MEEMER02_WA



30/08/2024 08:09:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

EMERGENCY BAG LIGHT (22L / 44X25X27)

EMERGENCY'S™

Ref. MEEMER02_WA

Practical and lightweight first aid bag in waterproof version.

Description

- It also has three outer pockets; the one at the front is divided with rubberized elastic bands
- High strength zip closures
- Reflective strips
- Washable case
- Reinforced handles and shoulder straps
- Dimensions 44 x 25 x 27 cm, capacity 22 L
- Material: Tarpaulin (laminated polyethylene)
- Weight 1.29 kg
- Sold empty

Features

CE Marking	NONE	Medical device	No
Bar code	3661040199191	HS Code	420212
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	25
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	27
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	44
Packed unit / Weight	1.3	Bare Unit / Weight (kg)	1.29
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 08:09:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr