

PLASTIC TROLLEY, 2 TRAYS 48,5 X 36,5 CM

Ref. OM483601



 29/08/2024 21:15:23 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PLASTIC TROLLEY, 2 TRAYS 48,5 X 36,5 CM

Ref. OM483601

Lightweight, handy and very robust pedestal table

Description

The quality of the casters and ergonomic handles ensure maximum comfort when travelling. - External dimensions: 48.5 x 36.5 x

H 72 cm, including handles

-Dimensions of the trays: 48 x 36.5 cm

-Maximum load: 100 kg per trolley and 60 kg per level

-Color: Light gray

-Materials: Aluminum pillar, PP top

- 4 PVC swivel casters, Ø 75 mm including 2 brake

-Weight: 4 kg

Features

CE Marking	NONE	Medical device	No
Bar code	3760206421296	HS Code	8716800000
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	41	Bare Unit / Width (cm)	36.5
Packed unit / Height (cm)	15	Bare Unit / Height (cm)	72
Packed unit / Length (cm)	56	Bare Unit / Length (cm)	48.5
Packed unit / Weight	5.2	Bare Unit / Weight (kg)	4
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	35
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 29/08/2024 21:15:23 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr