

## MALE CONNECTOR (FOR VACUUM PORT) FOR ASKIR 20 PACK OF 5'S



Ref. BCASPI164



25/08/2024 20:38:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## MALE CONNECTOR (FOR VACUUM PORT) FOR ASKIR 20 PACK OF 5'S



Ref. BCASPI164

### Features

CE Marking	CE	Medical device	No
Bar code	3661040160375	HS Code	90330090
Warranty (months)	12	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.01
Bare Unit / Weight (kg)	0.01		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

25/08/2024 20:38:40 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)