

SCHAEDEL TOWEL FORCEPS 9 CM

Ref. IP47009



27/08/2024 04:41:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SCHAEDEL TOWEL FORCEPS 9 CM



Ref. IP47009

Schaedel towel forceps 9 cm

Description

Schaedel towel forceps 9 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040026749	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.9
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	9
Packed unit / Weight	0.016	Bare Unit / Weight (kg)	0.016
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 04:41:48 - For information purposes only, can be subject to change