

DIAPHRAGME/GREY RETANNING RING FOR PULSE & DUAL PUSLE



Ref. 507937



30/08/2024 18:47:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DIAPHRAGME/GREY RETANNING RING FOR PULSE & DUAL PUSLE



Ref. 507937

Description

Spare part for Pulse® & Dual Pulse® stethoscopes

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446009687	HS Code	90181990
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.005
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 18:47:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr