

SET OF 5 SUNTECH MONOBLOC REUSABLE ARMBANDS ADULT PLUS (28-40CM) - BAYONET CONNECTOR

SUNTECH

Ref. DSU231016A



30/08/2024 09:43:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SET OF 5 SUNTECH MONOBLOC REUSABLE ARMBANDS ADULT PLUS (28-40CM) - BAYONET CONNECTOR

SUNTECH

Ref. **DSU231016A**

Description

Set of 5 Suntech adult plus nylon armbands (28-40cm)

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446013233	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	23	Packed unit / Height (cm)	9
Packed unit / Length (cm)	25	Packed unit / Weight	0.584
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


30/08/2024 09:43:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr