

SET OF 5 SUNTECH MONOBLOC DOUBLE TUBE REUSABLE ARMBANDS LARGE ADULT (31-40CM) - SCREW CONNECTOR

Ref. DSU231037

SUNTECH



30/08/2024 09:44:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

SET OF 5 SUNTECH MONOBLOC DOUBLE TUBE REUSABLE ARMBANDS LARGE ADULT (31-40CM) - SCREW CONNECTOR

SUNTECH

Ref. **DSU231037**

Description

Set of 5 Suntech burgundy nylon armbands.
Double tubing and screw connection, adult wide size 31-40cm

Features

CE Marking	CE	Medical device	Class I
Bar code	0120840935104114	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 09:44:38 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr