

NYLON ARMBAND SUNTECH ADUL (ROYAL BLUE) 23-33CM 2 TUBES WITHOUT FITTING

Ref. DSU231045

SUNTECH



 30/08/2024 09:42:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

NYLON ARMBAND SUNTECH ADUL (ROYAL BLUE) 23-33CM 2 TUBES WITHOUT FITTING

SUNTECH

Ref. **DSU231045**

Description

Suntech royal blue nylon armband individually.
Double tubing, without fitting, adult size 23-33 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446013967	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	14	Packed unit / Height (cm)	9
Packed unit / Length (cm)	19	Packed unit / Weight	0.448
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


30/08/2024 09:42:22 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr