

SET OF 20 SUNTECH MONOBLOC ARMBANDS SINGLE PATIENT LONG CHILD (12-19CM) - BAYONET CONNECTOR

Ref. DSU231062

SUNTECH



 30/08/2024 09:50:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SET OF 20 SUNTECH MONOBLOC
ARMBANDS SINGLE PATIENT LONG
CHILD (12-19CM) - BAYONET
CONNECTOR

SUNTECH

Ref. DSU231062

Description

Set of 20 Suntech "SOFT" single-patient armbands, green.
Child size long 12-19cm, rectus bayonnette connection.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446014100	HS Code	9018905040
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.857
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


30/08/2024 09:50:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr