

F3000 DISPOSABLE TEMP PROBE COVERS (25 BOXES/TRAY, 20 COVERS/BOX)

SUNTECH

Ref. DSU551467

CE
0413



 18/11/2024 21:15:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

F3000 DISPOSABLE TEMP PROBE COVERS (25 BOXES/TRAY, 20 COVERS/BOX)

SUNTECH

Ref. DSU551467

Description

25 boxes of 20

Features

CE Marking	CE0413	Medical device	Class I
Bar code	1088452116202	HS Code	9018.19.5500
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Packed unit / Height (cm)	10
Packed unit / Length (cm)	29	Packed unit / Weight	0.622
SPECIFIC STORAGE CONDITIONS			
OTHER FEATURES			
Packed unit / Quantity	1		

 18/11/2024 21:15:45 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr