

REPLACEMENT GUN FOR TEMPERATURE MODULE CT40

SUNTECH

Ref. DSU551470



 30/08/2024 09:53:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

REPLACEMENT GUN FOR TEMPERATURE MODULE CT40

SUNTECH

Ref. DSU551470

Features

CE Marking	CE0413	Medical device	Class IIa
Bar code	0110840935103813	HS Code	9018199560
Warranty (months)	3	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Packed unit / Height (cm)	8
Packed unit / Length (cm)	6	Packed unit / Weight	0.1
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 09:53:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr