

AUTOCLAVABLE CUFF ADULT EXTRA LONG LEG 102X13CM(BLADDER 84X9CM)

Ref. GHP552808

HpmTM



 30/08/2024 09:57:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr

AUTOCLAVABLE CUFF ADULT EXTRA LONG LEG 102X13CM(BLADDER 84X9CM)

HpmTM

Ref. GHP552808

Operating room equipment: double pocket autoclavable cuff for two-way electronic pneumatic tourniquet. Sold individually.

Description

Nylon cuff with single silicone bladder suitable for autoclaving. With male Luer lock connector and cap for cleaning and sterilisation.

Features

CE Marking	CE	Medical device	No
Bar code	3700446033729	HS Code	9018 90 2090
Warranty (months)	12	Country	Sweden
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.63
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 09:57:31 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr