

DISPOSABLE SINGLE CONICAL CUFF 76X10CM YELLOW (CARTON OF 10 PIECES)

Ref. GHP552901

Hpm™

CE
0123

DISPOSITIF MEDICAL
CLASSE
Is
MEDICAL DEVICE

GARANTIE
1
YEAR
GUARANTEE



📅 30/08/2024 09:59:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

DISPOSABLE SINGLE CONICAL CUFF 76X10CM YELLOW (CARTON OF 10 PIECES)



Ref. **GHP552901**

HPM armbands/bib shorts give you better access to the operating area thanks to a 2 m long TPU hose. The cuffs are made of high-quality material, which makes them soft and easy to apply, without the risk of latex allergy. Wide range of sizes with color coded for easy identification.

Description

ETO Sterilized. 2 m long hose made of TPU. Straight and Conical. No Latex. Wide range of sizes. Color coded. Customized protection sleeve included. CPC male as standard. Soft and easy to apply

Features

CE Marking	CE0123	Medical device	Class Is (sterile function)
Bar code	7350054600347	HS Code	9018 90 2090
Warranty (months)	12	Country	Sweden
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	27	Packed unit / Height (cm)	20
Packed unit / Length (cm)	64	Packed unit / Weight	3.174
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 09:59:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr