

SET OF 5 SUNTECH MONOBLOC INFANT REUSABLE ARMBANDS (8- 13CM) - SCREW CONNECTOR

Ref. DSU231000

SUNTECH



30/08/2024 09:35:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

SET OF 5 SUNTECH MONOBLOC INFANT REUSABLE ARMBANDS (8- 13CM) - SCREW CONNECTOR

SUNTECH

Ref. DSU231000

Description

Reusable nylon cuff sold by 5. Single tube, screw-on coupling.

Features

CE Marking	CE	Medical device	Class I
Bar code	0120840935100529	HS Code	9018905040
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.203
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

30/08/2024 09:35:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr