


## CLAMP FOR INFUSION POLE EGP 800

Ref. GHP552626

Hpm<sup>TM</sup>



 30/08/2024 10:58:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## CLAMP FOR INFUSION POLE EGP 800

Hpm™

Ref. GHP552626

## Features

|                             |               |                           |              |
|-----------------------------|---------------|---------------------------|--------------|
| CE Marking                  | NONE          | Medical device            | No           |
| Bar code                    | 3661040212449 | HS Code                   | 9018 90 2090 |
| Country                     | Sweden        |                           |              |
| DIMENSIONS AND WEIGHT       |               |                           |              |
| Packed unit / Width (cm)    | 11            | Packed unit / Height (cm) | nc           |
| Packed unit / Length (cm)   | 12            | Packed unit / Weight      | 0.242        |
| SPECIFIC STORAGE CONDITIONS |               |                           |              |
| Minimum temperature         | NA            | Max temperature           | NA           |
| Minimum humidity            | NA            | Maximum humidity          | NA           |
| OTHER FEATURES              |               |                           |              |
| Packed unit / Quantity      | 1             |                           |              |

 30/08/2024 10:58:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)