


MAGISTER STETHOSCOPE FOR CHILDREN / BLACK



Ref. 507217



30/08/2024 18:23:00 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

MAGISTER STETHOSCOPE FOR CHILDREN / BLACK



Ref. 507217

=====

Description

For pediatric auscultation. Stainless steel dual head chestpiece for pediatric use. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Equipped with a non-chill ring for little ones' comfort.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446009144	HS Code	90181990
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.375
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Pink	Packed unit / Quantity	1

30/08/2024 18:23:01 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr