

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 152 AND CLAMP/ADULT (M)



Ref. 337852



30/08/2024 18:04:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 152 AND CLAMP/ADULT (M)



Ref. 337852

Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Nihon Kohden patient monitors

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446000424	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Packed unit / Height (cm)	2
Packed unit / Length (cm)	21	Packed unit / Weight	0.128
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 30/08/2024 18:04:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr