

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 162 AND CLAMP/ADULT (M)



Ref. 337862



30/08/2024 18:05:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 162 AND CLAMP/ADULT (M)



Ref. 337862

Description

Single tube Clinic® cuff for Lian® + clamp compatible w/ Datascope, Drager, HP/ Agilent/ Philips, Marquette, Hellige, Nellcor, Siemens, Spacelabs patient

Features

CE Marking	NONE	Medical device	Class I
Bar code	3700446001940	HS Code	9018901000
Warranty (months)	6	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.13
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		


30/08/2024 18:05:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

 +33(0)4 42 90 31 31
 info@gsh-med.fr