


PLASTIC DISPOSABLE NOSE CLIP SIBELMED (X100)

Ref. DSSPIR03



 31/08/2024 01:07:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PLASTIC DISPOSABLE NOSE CLIP SIBELMED (X100)

Ref. DSSPIR03

Features

CE Marking	CE	Medical device	Class I
Bar code	8699089770403	HS Code	90192090
Warranty (months)	12	Country	Turkey
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	375
Bare Unit / Weight (kg)	310		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	50
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			

31/08/2024 01:07:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr