

SEAL FOR JAR 5000ML NEW HOSPIVAC



Ref. BCASPI123



 13/12/2024 22:10:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SEAL FOR JAR 5000ML NEW HOSPIVAC



Ref. BCASPI123

Features

CE Marking	NONE	Medical device	No
Bar code	0000000000741	HS Code	4202121900
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	17	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	17	Packed unit / Weight	0.024
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

13/12/2024 22:10:28 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr