

SUCTION UNIT CA-MI NEW HOSPIVAC 400 FLAC / SACH UU FLOVAC 2X2L FULL



Ref. BCASPI65



 13/12/2024 22:05:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 400 FLAC / SACH UU FLOVAC 2X2L FULL



Ref. BCASPI65

Description

FS 2 version with foot control and communicating jars

Designed to drain the patient's body fluids, gases, tissues and bones during or after a procedure. The state-of-the-art oil-free pump system guarantees excellent suction capabilities and maximum vacuum pressure in seconds, without any need for maintenance.

- The precise dashboard and the full range of accessories that come with it make it the ideal device for surgical aspiration and liposuction.
- Available in 3 versions: Basic, FS or full• Accessories for the Basic version:- Fluid collection jar with an anti-reflux valve

8x14 Silicone hoses (autoclavable) and conical connectors o 10-11-12mm

Antibacterial & hydrophobic filter (single patient)

Power cord with Schuko

plug• Accessories for the FS version:

Fluid collection jar with an anti-reflux

valve- 8x14 Silicone hoses (autoclavable) and conical connectors o 10-11-12mm- Antibacterial & hydrophobic filter (single patient)

Power cord with Schuko plug

Pedal

• Accessories for the Full version:

Fluid collection jar with an anti-reflux valve

8x14 Silicone hoses (autoclavable) and conical connectors o 10-11-12mm

- Antibacterial & hydrophobic filter (single patient)

- Power cord with Schuko plug,

Pedal

Electronic Transfer System between collection jars

• Technical characteristics

:- Engine: Oil-free piston pump requiring no maintenance

Power supply: 230V-50Hz

13/12/2024 22:05:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SUCTION UNIT CA-MI NEW HOSPIVAC 400 FLAC / SACH UU FLOVAC 2X2L FULL



Ref. BCASPI65

Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040139517	HS Code	90189084
Warranty (months)	24	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	20
Bare Unit / Weight (kg)	20		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

13/12/2024 22:05:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr