

LUER BONE RONGEUR STRAIGHT 17 CM 8 MM

Ref. IP48017



10/01/2025 02:15:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

LUER BONE RONGEUR STRAIGHT 17 CM 8 MM



Ref. IP48017

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040026992	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	17
Packed unit / Weight	0.147	Bare Unit / Weight (kg)	0.147
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

10/01/2025 02:15:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr