

Réf. BCASPI17



CE
0123

DISPOSITIF MEDICAL
MEDICAL DEVICE

CLASSE
IIa

Réf. BCASPI17

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3661040139036
Guarantee (month):	36	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :	4,400	Weight packaging (kg) :	4,400
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/10/2022 12:39 - This article is available on our website : www.gsh-med.fr with the following reference BCASPI17

All characteristics are given for information only and may be changed