



MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507209



09/01/2025 23:21:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

MAGISTER STETHOSCOPE FOR ADULT / BLACK



Ref. 507209

=====

Description

For general-purpose auscultation. Stainless steel dual head chestpiece for adults. Multifrequency diaphragm. Angled binaural with built-in spring and flexible, screw-on tips. Ideal grip and high performance make it the perfect product for general and specialist use.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446004002	HS Code	9018.90.20.00-9
Warranty (months)	60		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.398
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Colors	Green	Packed unit / Quantity	1

09/01/2025 23:21:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr