

BARRAQUER 'COLIBRI' TYING FORCEPS - 0.5 MM

Ref. IP52101



09/01/2025 21:48:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

BARRAQUER 'COLIBRI' TYING FORCEPS - 0.5 MM



Ref. IP52101

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040027791	HS Code	90189090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	7.6
Packed unit / Weight	0.008	Bare Unit / Weight (kg)	0.008
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

09/01/2025 21:48:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr