


## STAINLESS STEEL EXAMINATION COUCH (QUICK ASSEMBLING) BLACK



Ref. OSA003\_N



14/01/2025 22:53:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## STAINLESS STEEL EXAMINATION COUCH (QUICK ASSEMBLING) BLACK



Ref. OSA003\_N

### Quick Mount System

### Description

- Rapid assembly system
- Square stainless steel tube frame
- Adjustable headboard with automatic high resistance rack
- Maximum weight: 150 kg
- Headboard and seat covered with 5 cm thick very high density foam (30 kg/m3), covered by fire-resistant washable plasticised fabric (fire resistance class: M1)
- Removable roll holder (for 50 cm wide rolls) that can be located at head or feet

### OPTIONS

Ref. OTA0047 Heel support without clamps (pair)

Ref. OTA0041 Leg holders without clamps (pair)

Ref. OTA0044 Fixing clamps (pair)

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040165844	HS Code	9402900000
Warranty (months)	36	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Bare Unit / Width (cm)	63
Packed unit / Height (cm)	68	Bare Unit / Height (cm)	70
Packed unit / Length (cm)	190	Bare Unit / Length (cm)	190
Packed unit / Weight (kg)	27.5	Bare Unit / Weight (kg)	24.5

For information purposes only, can be subject to change

### SPECIFIC STORAGE CONDITIONS

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)