

KIT EXAMINATION COUCH+STOOL+ MOUNTING STEP 1 STEP S/S-BLACK - HOLTEx



Ref. OSKIT001



16/10/2024 21:18:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

KIT EXAMINATION COUCH+STOOL+ MOUNTING STEP 1 STEP S/S-BLACK - HOLTEx



Ref. OSKIT001

Complete stainless steel cabinet

Description

Delivered in one box. Set of furniture, consisting of:- 1 stainless steel sofa black color (ref OSA003_N)- 1 stool black color (ref OS007_N)- 1 step (ref OS0051)

The sofa with quick assembly system:

- Structure in square stainless steel tube
- Adjustable headrest with automatic rack high resistance
- Head and seat padded with very high density foam (30 kg / m3), whose thickness is 5 cm, covered with washable plasticized fabric, fireproof (fire resistance classification: M1)
- Removable roller holder (for rolls 50 cm wide) positionable at the head or feet
- Max load: 150 kg

The stool with 5 swivel casters,

- structure 5 branches ABS, 58cm in diameter- 1 gas piston for a variable height 53 to 70 cm- 1 foam seat 30 kg / M3, diameter 43 cm, with screwed steel plate
- Plasticized coating, washable, fireproof classification M1- Max load: 150 kg

The step 1 step

- stainless steel and non-slip plastic coating
- Load: 150 kg

Stainless steel or Epoxy: What's the difference?

Epoxy: the economic choice. The steel is simply coated with a resin to delay corrosion.

Stainless steel: it is steel modified to become completely STAINLESS steel, it does not rust.

And did you know?

Today, more than 90% of stainless steel is made from recycled steel, it's good for the environment!

📅 16/10/2024 21:18:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr

KIT EXAMINATION COUCH+STOOL+ MOUNTING STEP 1 STEP S/S-BLACK - HOLTEx



Ref. OSKIT001

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040168371	HS Code	9402900000
Warranty (months)	24	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Packed unit / Height (cm)	68
Packed unit / Length (cm)	190	Packed unit / Weight	36.7
Bare Unit / Weight (kg)	33.7		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

16/10/2024 21:18:16 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr