

TÜBINGEN SUTURE FORCEPS 10CM STRAIGHT 1X2 TEETH

Ref. IP52913



27/08/2024 05:10:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

TÜBINGEN SUTURE FORCEPS 10CM STRAIGHT 1X2 TEETH



Ref. IP52913

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040027968	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	0.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	10
Packed unit / Weight	0.014	Bare Unit / Weight (kg)	0.014
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 05:10:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr