

## SCHWEIGGER CAPSULAR FORCEPS MULTI TEETH ANGULATED

Ref. IP53110



27/08/2024 05:13:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## SCHWEIGGER CAPSULAR FORCEPS MULTI TEETH ANGULATED



Ref. IP53110

### Features

|                             |               |                         |           |
|-----------------------------|---------------|-------------------------|-----------|
| CE Marking                  | CE            | Medical device          | Class Ir  |
| Bar code                    | 3661040028071 | HS Code                 | 9018.9090 |
| Warranty (months)           | 24            | Country                 | Pakistan  |
| DIMENSIONS AND WEIGHT       |               |                         |           |
| Packed unit / Width (cm)    | nc            | Bare Unit / Width (cm)  | 4         |
| Packed unit / Height (cm)   | nc            | Bare Unit / Height (cm) | 1.8       |
| Packed unit / Length (cm)   | nc            | Bare Unit / Length (cm) | 13.1      |
| Packed unit / Weight        | 0.032         | Bare Unit / Weight (kg) | 0.03      |
| SPECIFIC STORAGE CONDITIONS |               |                         |           |
| Minimum temperature         | NA            | Max temperature         | NA        |
| Minimum humidity            | NA            | Maximum humidity        | NA        |
| OTHER FEATURES              |               |                         |           |
| Packed unit / Quantity      | 1             |                         |           |

27/08/2024 05:13:05 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France  
[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31  
[info@gsh-med.fr](mailto:info@gsh-med.fr)