

KOERTE RETRACTOR 26 X 40 MM - BLADE 235 MM BLUNT

Ref. IE11600



27/08/2024 00:53:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

KOERTE RETRACTOR 26 X 40 MM - BLADE 235 MM BLUNT



Ref. IE11600

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040084725	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	4.3
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.9
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	24.6
Packed unit / Weight	0.11	Bare Unit / Weight (kg)	0.107
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 27/08/2024 00:53:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr