

CARTRIGE SYRINGES WITHOUT ASPIRATION LATERAL LOADING 1.8 CC



Ref. IDSECSL1



27/08/2024 00:32:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

CARTRIGE SYRINGES WITHOUT ASPIRATION LATERAL LOADING 1.8 CC



Ref. IDSECSL1

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040015958	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	11.3
Packed unit / Weight	0.091	Bare Unit / Weight (kg)	0.091
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 00:32:54 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr