

KNAPP RETRACTOR, 14CM

Ref. IE41100



27/08/2024 01:27:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

KNAPP RETRACTOR, 14CM



Ref. IE41100

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040017709	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	4.4
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	2.1
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	17.2
Packed unit / Weight	0.021	Bare Unit / Weight (kg)	0.021
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 01:27:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr