

JOHN HOPKINS AORTIC CLAMP

Ref. IC30200



26/08/2024 22:04:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr

JOHN HOPKINS AORTIC CLAMP



Ref. IC30200

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040010946	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	5	Bare Unit / Width (cm)	5.5
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2.7
Packed unit / Length (cm)	25	Bare Unit / Length (cm)	25.5
Packed unit / Weight	0.122	Bare Unit / Weight (kg)	0.122
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 22:04:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr