

RIB SHEAR 20 CM DOYEN TYPE

Ref. IC32420



10/01/2025 04:16:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
 www.gsh-med.fr

+33(0)4 42 90 31 31
 info@gsh-med.fr

RIB SHEAR 20 CM DOYEN TYPE



Ref. IC32420

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040011028	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	8.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.4
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	21
Packed unit / Weight	0.253	Bare Unit / Weight (kg)	0.253
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

10/01/2025 04:16:47 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France
www.gsh-med.fr

+33(0)4 42 90 31 31
info@gsh-med.fr